



## ArieSpike Bag

Perforatori per il travaso  
e la ricostituzione dei farmaci in sacche

Non vented spike for withdrawal  
and reconstitution of drugs in bags



## Protezione degli operatori sanitari con dispositivi non ventilati privi di ago per il travaso da sacche multidose

- Sicuro. Il sistema di aggancio brevettato, ideato per il fissaggio delle alette alle sacche BAXTER "Viaflo", consente di evitare sconnessioni accidentali. Il perforatore, disegnato per adattarsi perfettamente al canale di connessione alla sacca, garantisce una tenuta sicura nonché il completo recupero della soluzione senza lasciare residui.
- Da utilizzare in ambiente a contaminazione controllata durante la fase di prelievo del farmaco da sacche multidose preconfezionate.
- Facilità di utilizzo. Sistema senz'ago monolume per il controllo delle preparazioni sterili.
- Il tappo luer lock o la valvola autosigillante permettono lo stoccaggio sicuro del preparato consentendo il prelievo sequenziale da sacche multidose.
- La capsula sulla punta del perforatore protegge gli operatori durante la fase di manipolazione del dispositivo.

## Protects healthcare professionals with needle-less non vented access for multi-dose bags withdrawals

- Safe. The patented connection system, studied to be connected to the "wings" of Baxter "Viaflo" bags, avoids accidental disconnection of the spike from the bags. The spike, studied to perfectly adapt to the bag withdrawal channel of the bags, allows safety connection and complete recovery of the drug from the bags.
- Use in pharmacy "cleanroom" or nursing patient care environment during solution withdrawal from prepackaged multi-dose bags.
- Easy to use. Needle-less single lumen system for quality control of sterile drug preparation.
- Luer lock cap or needle less valve permits safe storage and allows sequential withdrawal from multi-dose bags.
- Capped spike protects health care workers during handling.

Codice Ref. Nr.	Filtro aria Air Vent	Filtro soluzione Solution Filter	Perforatore Spike	Valvola autosigillante Needle less valve	Senza PVC/Lattice PVC free/Latex free	Materiale Material
CH010021	No		ABS		•	ABS
BS051004	No		ABS	•	•	ABS
BS051009	No	5µ	ABS	•	•	ABS