



 UV-Barrier

Sacche trasparenti per la protezione
delle sostanze fotosensibili

Photoprotective bags for protection
of light sensitive solutions



Sacche per la protezione dei contenitori (sacche, flaconi, siringhe) contenenti sostanze fotosensibili.

- Materiale di fabbricazione: Polietilene LD. Priva di lattice e PVC
- Prodotte con speciale pigmento ad effetto barriera per fonte luminosa con lunghezza d'onda compresa tra 200 e 600 nm.
- Il filtro ottico consente di bloccare selettivamente le lunghezze d'onda responsabili della inattivazione delle sostanze fotosensibili, mantenendo al tempo stesso la trasparenza della sacca.
- Queste sacche sono in grado di proteggere le sostanze fotosensibili da quei fenomeni di fotodegradazione selettiva che porta alla diminuzione se non in alcuni casi alla scomparsa degli effetti terapeutici.
- Ogni sacca della gamma UV-Barrier è corredata di striscia adesiva per la chiusura, foro di aggancio, asola di uscita per la connessione al set di somministrazione.
- Classe di appartenenza: 1 NON STERILE
- Codice CND: A0899.

Performance

Accessorio per la preparazione e protezione dei farmaci.

Protective bags for vial, bags and syringes containing light sensitive solutions.

- Materials: polyethylene LD. Latex free – pvc free
- Produced with a special barrier effect for light source with wave length between 200 and 600 nm.
- The optic filter allows to block selectively the wave length responsible for the inactivation of light sensitive solutions. The transparency of the material allows to check the liquid inside.
- These bags are studied to protect the light sensitive solutions from those problems of selective photo degradation which could lead to a partial or total disappearance of therapeutic effects.
- Each bag of the UV-Barrier range is produced with a special adhesive sealing stripe, an open site for the administration sets and a hole to hang them on a pole.
- Class: I NON STERILE
- CND Ref N°: A0899

Performance

Liquid solutions protection device



Codice Ref. Nr.	Misure Sizes	Indicazione Use for	Trasparente Transparent	Senza PVC PVC free	Chiusura adesiva Adhesive sealing	Senza PVC/Lattice PVC free/Latex free	Materiale Material
CH000001	95 x 360 mm.	siringhe / syringes	•	•	•	•	PE
CH000002	100 x 230 mm.	flaconi/sacche 100 ml. vials/bags 100 ml.	•	•	•	•	PE
CH000003	130 x 300 mm.	flaconi/sacche 250 ml. vials/bags 250 ml	•	•	•	•	PE
CH000004	165 x 350 mm.	flaconi/sacche 500/1000 ml. vials/bags 500/1000 ml	•	•	•	•	PE
CH000007	350 x 440 mm.	sacche 3000/2000 ml. bags 3000/2000 ml	•	•	•	•	PE