

DEVICES FOR UROLOGICAL WASHING

ARIES S.r.l.
Biomedical Devices



ARIES S.r.l.
Biomedical Devices

DEVICES FOR UROLOGICAL WASHING

DRAWING	DESCRIPTION	CODE	PACK. QTY.
	BioDrip 1-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Drip chamber: in phthalate-free PVC, height 10 cm, volume 20 ml. ABS roller clamp. Tubing: Dehp free PVC Ø 4.8 x 6.8. Terminal fitting: stepped catheter cone adapter in Dehp free PVC equipped with silicone connection tube. Length 240 cm	UR010001.	50
	BioDrip 2-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Clamps: in polypropylene with snap closure and opening. Drip chamber: in Dehp free PVC, height 10 cm, volume 20 ml. ABS roller clamp. Tubing: Dehp free PVC Ø 4.8 x 6.8. Terminal fitting: stepped catheter cone adapter in Dehp free PVC equipped with silicone connection tube. Length 240 cm.	UR010002.	40



ARIES S.r.l.
Biomedical Devices

DEVICES FOR UROLOGICAL WASHING

DRAWING	DESCRIPTION	CODE	PACK. QTY.
	BioDrip 2-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Clamps: in polypropylene. Drip chamber: in Dehp free PVC, height 10 cm, volume 20 ml. Clamp: roller clamp in ABS. Tubing: Dehp Free PVC Ø 4.8 x 6.8. Terminal fitting: LL/ Male in ABS with protection capsule. Version with: • Length 220 cm. • Length 300 cm.	UR010006 UR010014	40
	BioDrip 2-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Clamps: in polypropylene. Drip chamber: in Dehp free PVC, height 10 cm, volume 20 ml. Clamp: roller clamp in ABS. Tubing: Dehp Free PVC Ø 4.8 x 6.8. Terminal fitting: LL/ Male in ABS with protection capsule. Length 220 cm.	UR010007	40


ARIES S.r.l.
Biomedical Devices

DEVICES FOR UROLOGICAL WASHING

DRAWING	DESCRIPTION	CODE	PACK. QTY.
	BioDrip 2-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Clamps: in polypropylene. Drip chamber: in Dehp free PVC, height 13 cm with 2-inlet lid, Ø30. Tubing: Dehp Free PVC Ø 4.8 x 6.8. Terminal fitting: PVC tube. Length 300 cm.	UR010012	30
	BioDrip 5-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Clamps: in polypropylene with snap closure and opening. Drip chamber: in PVC, height 10 cm, volume 20 ml. Clamp: roller clamp in ABS. Tubing: Dehp Free PVC Ø 4.8 x 6.8. Terminal fitting: stepped catheter cone adapter in Dehp free PVC equipped with connection tube. Length 240 cm.	UR010005.	40

ARIES S.r.l.
Biomedical Devices

DEVICES FOR UROLOGICAL WASHING

DRAWING	DESCRIPTION	CODE	PACK. QTY.
	BioDrip 8-way urological infusion set. Spike: 1-way in ABS without air intake for use with bags. Clamps: in polypropylene with snap closure and opening. Drip Chamber: Dehp Free PVC Ø 16. Clamp: roller clamp in ABS. Tubing: Dehp Free PVC Ø 4.8 x 6.8. Terminal fitting: stepped catheter cone adapter in Dehp free PVC equipped with connection tube. Length 275 cm.	UR010003	20
	IV stand Adjustable height from 1500 mm to 2300 mm using a foot-operated stainless steel hydraulic pump 5-star stainless steel base Moving castors Ø 80 mm of which two with brake 4 pigtail hooks Payload: 10 Kg each hook Weight: 23 Kg		Sold individually

ARIES S.r.l.
Biomedical Devices

NOTES:

ARIES S.r.l.
Biomedical Devices

NOTES: