

CONTAINMENT AND ADMINISTRATION SYSTEMS FOR PARENTERAL NUTRITION

ARIES S.r.l.
Biomedical Devices



ARIES S.r.l.
Biomedical Devices

CONTAINMENT AND ADMINISTRATION SYSTEMS FOR PARENTERAL NUTRITION

DRAWING	DESCRIPTION	CODE	PACK. QTY.
	Eval Parenteral nutrition systems. Bag: EVA knurled configuration to prevent the internal walls from collapsing. Spike: 1-way in ABS without air intake. Injection site: latex-free pierceable membrane. Pierceable fitting: in soft PVC. 3-way filling set: 2-way spike, slide clamp, length of each branch 30 cm. Tubing: Dehp free PVC length 90 cm. • Capacity 500 ml. • Capacity 1000 ml. • Capacity 2000 ml. • Capacity 3000 ml.	NP011003. NP011004. NP011005. NP011006.	30
	Eval Parenteral nutrition systems. Bag: EVA knurled configuration to prevent the internal walls from collapsing. Spike: 1-way in ABS without air intake. Injection site: latex-free pierceable membrane. Pierceable fitting: in soft PVC. 3-way filling set: 2-way spike, slide clamp, length of each branch 30 cm. Tubing: Dehp free PVC length 90 cm. • Capacity 500 ml. • Capacity 1000 ml. • Capacity 2000 ml. • Capacity 3000 ml.	NP011103 NP011102 NP011101 NP011100	40 30 30 30

ARIES S.r.l.
Biomedical Devices

CONTAINMENT AND ADMINISTRATION SYSTEMS FOR PARENTERAL NUTRITION

DRAWING	DESCRIPTION	CODE	PACK. QTY.
	Eval Parenteral nutrition systems. Bag: EVA knurled configuration to prevent the internal walls from collapsing. Injection site: latex-free pierceable membrane. Pierceable fitting: in soft PVC. Fitting: in rigid PVC complete with luer lock connection. • Capacity 250 ml. • Capacity 500 ml. • Capacity 1000 ml.	NP011096. NP011095. NP011094.	50
	Eval Parenteral nutrition systems. Bag filling system. 2-way spikes: in ABS with air intake and 3µm water-repellent antibacterial filter. Clamp: in Polyethylene Tubing: Transparent soft DEHP FREE PVC Ø 3 x 4.1 Length 70 cm. 3-way joined fitting in rigid PVC One-way valve: in polycarbonate with Luer Lock connection Needle: with drawn steel bar and Luer Lock fitting Syringe: 50/60 ml in transparent polypropylene	NP021000	40

ARIES S.r.l.
Biomedical Devices

NOTES:

NOTES: